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**Verification of qRT-PCR assay for detection of SARS-CoV-2**

*Version 2; 20-05-2021: concentration dPCR copies SARS-CoV-2 updated and some additional information added on quantitation limit. This only affects the description of the absolute minimum amount that can be detected. The test remains having the highest sensitivity possible in national and international comparisons.*

**Method:** qPCR – One step RT-PCR with Fast Virus Master Mix

**Manufacturer/Origin:** Biolegio

**Type of validation:** verification of improved primers/probe set specific for detection of SARS-CoV-2

**Application test:** patient diagnostics, detection of SARS-CoV-2 in clinical material

**Conclusion:** Based on the results, we conclude that using the improved primers/probe set the RT-PCR meets the predefined criteria of the following parameters: measurement trueness, measurement accuracy, quantitation limit, analytical sensitivity and detection limit. The RT-PCR is more sensitive than with the original primers/probe set and now has similar sensitivity as the E-gene RT-PCR. Therefore the improved primers/probe set has been accepted for implementation in routine diagnostics of SARS-CoV-2.

**Author:** 5.1.2e

**Date:** 6-7-2020

**Introduction & purpose**

At this moment when a request is made for corona diagnostics the sample is tested on 2 genes of SARS-CoV-2. The E-gene and RdRp-gene are detected in 2 separate PCRs. The E-gene PCR is more sensitive than the RdRp-gene PCR. With an adjustment of the primers and probe set for the RdRp-gene, the idea is that the sensitivity of RdRp-gene PCR improves, so that it provides additional confirmation whether a sample is SARS-CoV-2 positive. Especially when the viral load is low.

**Material & method**

Wobbles are removed from the original primers to match the target sites in SARS-CoV-2 exactly, the reverse primer is extended to better match the T<sub>m</sub> of the forward primer and the probe has been extended by 1 nt for increased T<sub>m</sub> and better binding compared to the original one for an annealing/extension temperature of 60°C (adjustment based on work from 5.1.2e JBZ, personal communication). The adjusted primers and probe were analysed against a representative (geographic origin and distributed over time in the pandemic) set of sequences for match; no mismatched nucleotides were found. Samples for testing were extracted on the MagNa Pure 96 instrument (Roche) using the MagNa Pure 96 DNA and Viral NA Small volume kit (Roche) and eluted in a volume of 50 µl. RdRp-assay primers and probes as described in the provided protocol by Corman et al.[1] (Table 1) were used to compare the performance of the new specific SARS-CoV2 RdRp gene primers and probe (Table 2). Reagents and temperature protocol are described in Tables 3 and 4 respectively.

**Table 1:** Original RdRp-gene primers and probes Corman et al.[1]

| Oligonucleotide | Sequence                           | Label    | Target     |
|-----------------|------------------------------------|----------|------------|
| RdRp_SARSr-F    | GTG ARA TGG TCA TGT GTG GCG G      |          | Sarbeco    |
| RdRp_SARSr-R    | CAR ATG TTA AAS ACA CTA TTA GCA TA |          | Sarbeco    |
| RdRp_SARSr-P1   | CCA GGT GGW ACR TCA TCM GGT GAT GC | TXR-BHQ2 | Sarbeco    |
| RdRp_SARSr-P2   | CAG GTG GAA CCT CAT CAG GAG ATG C  | FAM-BHQ1 | SARS-CoV-2 |

**Table 2:** Adjusted primers and probe specific for SARS-CoV2 RdRp gene

| Oligonucleotide | Sequence                                                     | Label    |
|-----------------|--------------------------------------------------------------|----------|
| RdRp_SARS-F2    | GTG <b>AAA</b> TGG TCA TGT GTG GCG G                         |          |
| RdRp_SARS-R2    | <b>CAA</b> ATG TTA <b>AAA</b> ACA CTA TTA GCA TAA <b>GCA</b> |          |
| RdRp_SARS-P2.2  | <b>CCA</b> GGT GGA ACC TCA TCA GGA GAT GC                    | FAM-BHQ1 |

**Table 3:** Composition reagents mixture for RdRp-gene SARS-CoV-2 RT-PCR. The primers and probes obtained from Biolegio were premixed at a final concentration of 10 µM primers and 5 µM probe.

|                           |       |
|---------------------------|-------|
| RdRp qRT PCR              | µl    |
| 4x Taqman Fast Virus MM   | 5     |
| Primers/probe Mix (10 µM) | 3     |
| PCR grade water           | 7     |
| Specimen nucleic acid     | 5     |
| Total volume              | 20 µl |

**Table 4:** Amplification temperature protocol.

| PCR Program               | Segment number | Temp Target (°C) | Hold Time (sec.) | Slope (°C/sec.) | Acquisition mode |        |
|---------------------------|----------------|------------------|------------------|-----------------|------------------|--------|
| Reverse Transcription     | 1              | 50               | 900              | EXTERNAL        |                  |        |
| Denaturation/Inactivation | 1              | 95               | 120              | EXTERNAL        |                  |        |
| Denaturation              | 1              | 95               | 60               | 4.4             | None             | LC 480 |
| Amplification             | 1              | 95               | 10               | 4.4             | None             |        |
| (cycles:50)               | 2              | 60               | 30               | 2.2             | Single           |        |
| Cooling                   | 1              | 40               | 30               | 4.4             | None             |        |

**Analytic parameters:**

The detection qPCR is not used as a pure quantitative method. Results are scored positive/negative to confirm if a patient is infected with SARS-CoV-2. Ct values are used as semi-quantitative information when useful for interpretation of the result. Only those parameters are tested that need verification compared to the original validation

|                                                                                                                                           |                                                           |
|-------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------|
| <u>Measurement uncertainty</u> (Meetonzekerheid):                                                                                         | – not tested                                              |
| <u>Measurement trueness</u> (Juistheid):                                                                                                  | – tested: Saliva panel (page 5)                           |
| <u>Measurement accuracy</u> (Accuraatheid):                                                                                               | – tested: Saliva panel (page 5)                           |
| <u>Measurement precision</u> including measurement repeatability (Precisie inclusief herhaalbaarheid):                                    | – not tested                                              |
| <u>Analytical sensitivity</u> (Analytische sensitiviteit):                                                                                | – tested: LOD95 determination (page 4)                    |
| <u>Analytical specificity</u> , including interfering substances (Analytische specificiteit met inbegrip van interfererende substanties): | – not tested                                              |
| <u>Detection limit</u> (Detectielimiet, bij kwantitatieve methoden):                                                                      | – tested: LOD95 determination (page 4)                    |
| <u>Quantitation limit</u> (Kwantificatielimiet, bij kwantitatieve methoden):                                                              | – tested: 10-fold dilution of quantified control (page 3) |
| <u>Measuring interval</u> (Meetinterval, kwantitatieve methoden):                                                                         | – not tested                                              |

**Clinical parameters:**

|                                                              |              |
|--------------------------------------------------------------|--------------|
| <u>Diagnostic sensitivity</u> (Diagnostische Sensitiviteit): | – not tested |
| <u>Diagnostic specificity</u> (Diagnostische Specificiteit): | – not tested |

## Results

### Quantitation limit

The SARS-CoV-2 RNA used for this verification is an 1:100 dilution of the hCoV-19/Netherlands/NoordBrabant\_10003/2020 (4732000117 P2) stock obtained from cell culture. To establish PCR efficiency we first ran a triplicate 10-fold dilution series of viral RNA for each assay. Viral RNA was isolated by using the MagNA Pure LC Total Nucleic Acid Isolation Kit (Roche). After Nucleic Acid extraction the concentration is  $5.58 \times 10^5$  digital copies RdRp-gene SARS-CoV-2 positive strand RNA/ $\mu$ l extract.

**Table 5:** 10-fold dilution series both primers and probe sets

| Dilution         | Current nCoV detection |          |          | New nCoV detection |          |          |
|------------------|------------------------|----------|----------|--------------------|----------|----------|
|                  | I                      | II       | III      | I                  | II       | III      |
| SARS-CoV-2 10-1* | 23.81                  | 24.31    | 24.55    | 22.07              | 22.17    | 22.04    |
| SARS-CoV-2 10-2  | 27.93                  | 28.16    | 28.27    | 25.77              | 26.04    | 25.94    |
| SARS-CoV-2 10-3  | 30.63                  | 31.17    | 30.9     | 28.97              | 28.95    | 28.89    |
| SARS-CoV-2 10-4  | 33.23                  | 33.27    | 32.9     | 32.41              | 32.45    | 32.5     |
| SARS-CoV-2 10-5  | 34.91                  | 34.92    | 34.55    | 34.92              | 36.2     | 34.97    |
| SARS-CoV-2 10-6  | Negative               | Negative | Negative | Negative           | 37.52    | Negative |
| SARS-CoV-2 10-7  | Negative               | Negative | Negative | 37.25              | 37.19    | Negative |
| SARS-CoV-2 10-8  | Negative               | Negative | Negative | Negative           | Negative | Negative |

\* Starting concentration is  $5.58 \times 10^4$  RNA copies/ $\mu$ l extract.

We determined the slope in Excel and defined the required levels for PCR efficiency (E) and R<sup>2</sup> as >95% and >0.95, respectively. The new specific SARS-CoV-2 RdRp gene primers and probe meet the requirements (Figure 1, Table 6).


**Figure 1:** Linearity and amplification efficiency both primers and probes sets.

**Table 6:** The amplification efficiency and linearity calculated using the slope of the targets in figure 1 (calculated in Excel)

|              | Slope  | $E = 10^{-1/\text{slope}}$ | $E(\%) = 100 \times [10^{-1/\text{slope} - 1}]$ | R <sup>2</sup> |
|--------------|--------|----------------------------|-------------------------------------------------|----------------|
| Theorie      | -3.322 | 2.00                       | 100                                             | 0.9500         |
| Current nCoV | -2.615 | 2.41                       | 141.19                                          | 0.9739         |
| New nCov     | -3.308 | 2.01                       | 100.60                                          | 0.9981         |

### Analytical sensitivity and Detection Limit

Next, we ran four replicates of a 2-fold dilution series (diluted in yeast carrier RNA in water) to determine the LOD95% by Probit analysis using SPSS Statistics (IBM, version 24) (Figures 2 and 3).


**Figure 2:** The results in percentages from the four replicates of the 2-fold dilution series against copies RdRp RNA/μl extract.

### LOD95% with upper and lower bound


**Figure 3:** The LOD95% in RdRp RNA copies/μl extract of the current nCoV-RdRp-gene and the new nCoV-RdRp-gene RT-qPCR.


### Measurement trueness and Measurement accuracy

The measurement trueness and accuracy has been tested using the saliva EQA panel. The panel specimens were extracted in triplicate before and after storage at -80 °C. The sensitivity of the new nCoV PCR is better than the current test. In the current test the 10<sup>-7</sup> dilution, virus particles equivalent to 1395 RdRp-gene digital PCR copies/ml saliva, is negative. In the new nCoV test the 10<sup>-7</sup> dilution is positive (Table 7 and Figure 4). A freeze/thaw cycle did not result in reduction of detectable RNA.

**Table 7:** Saliva panel with cp-value\*

| Name                 | Dilution Stockvirus<br>4732000117 P2 | Cp nCoV new |       |       | Cp nCoV current |       |       | Comment               |
|----------------------|--------------------------------------|-------------|-------|-------|-----------------|-------|-------|-----------------------|
| Sen.Saliva CoV20-1   | 10-4                                 | 25.46       | 25.48 | 25.21 | 25.86           | 26.02 | 25.75 | Before -80 °C storage |
| Sen.Saliva CoV20-2   | 10-7                                 | 32.01       | 31.91 | 32.04 |                 |       |       | Before -80 °C storage |
| Sen.Saliva CoV20-3   | 10-5                                 | 28.84       | 28.7  | 28.62 | 27.76           | 27.62 | 28.44 | Before -80 °C storage |
| Sen.Saliva CoV20-4   | negative                             |             |       |       |                 |       |       | Before -80 °C storage |
| Sen.Saliva CoV20-5   | 10-6                                 | 29.42       | 29.15 | 30.75 | 28.19           | 28.02 | 28.72 | Before -80 °C storage |
| Sen.Saliva CoV20-6   | 10-8                                 |             |       | 32.55 |                 |       |       | Before -80 °C storage |
| Sen.Saliva CoV20-1.1 | 10-4                                 | 25.2        | 25.94 | 26    | 25.86           | 26.21 | 26.23 | After -80 °C storage  |
| Sen.Saliva CoV20-2.1 | 10-7                                 | 31.82       | 31.98 | 31.95 |                 |       |       | After -80 °C storage  |
| Sen.Saliva CoV20-3.1 | 10-5                                 | 29.7        | 28.18 | 28.82 | 28.11           | 27.56 | 28    | After -80 °C storage  |
| Sen.Saliva CoV20-4.1 | negative                             |             |       |       |                 |       |       | After -80 °C storage  |
| Sen.Saliva CoV20-5.1 | 10-6                                 | 30.98       | 30.88 | 31    | 29.46           | 29.12 | 28.94 | After -80 °C storage  |
| Sen.Saliva CoV20-6.1 | 10-8                                 |             |       |       |                 |       |       | After -80 °C storage  |

\* Concentration of the dilution series starts at a virus concentration equivalent to  $1.40 \times 10^6$  (10-4 dilution) dPCR RdRp-gene copies/ml in specimen.

**Figure 4:** E-gene and RdRp-gene qRT-PCR results of saliva panel before and after -80 °C storage

### Conclusion

Based on the results of the measurement trueness, measurement accuracy, quantitation limit, analytical sensitivity and detection limit the test can be implemented to detect SARS-CoV-2 in clinical specimens.

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**Validation of qRT-PCR assays for detection of novel coronavirus 2019-nCoV**

*Version 2; 20-05-2021: number of copies SARS-CoV-1/μl harmonized and some additional information added on quantitation limit. This only affects the description of the absolute minimum amount that can be detected. The test remains having the highest sensitivity possible in national and international comparisons.*

**Methode:** qPCR – Onestep RT-PCR met behulp van Fast Virus Master Mix

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**Auteur:** 5.1.2e

**Datum:** 22-1-2020

**Introduction & goal****Background surveillance acute respiratory infections**

The clinical specimens for the clinical validation of the E-assay and GdRp-assay were selected from the biobank containing combined nose swab and throat swab specimens in virus transport medium from the national surveillance of influenza-like illness (ILI) and acute respiratory infections (ARI) in the Netherlands. General practitioners from the Nivel Netherlands institute for health services research Primary Care Database report the incidence of ILI on a weekly basis and take nasopharyngeal and oropharyngeal swabs from patients consulting with ILI or ARI. These specimens are submitted to the National Institute for Public Health and the Environment (RIVM) for in house developed qRT-PCR diagnosis of influenza virus, respiratory syncytial virus, rhinovirus and enterovirus. When questions arise on unsuspected increases in ILI incidence accompanied by a low detection rate of pathogens the FTD Respiratory pathogens 21 assay (Fast Track Diagnostics) is used.

**Material & methods****Selection of specimens**

One-hundred specimens have been selected from the surveillance biobank (Table 1). For influenza viruses, RSV and rhinovirus about ten specimens with Ct values below 30 and representative for the range of Ct values in clinical practice. For enterovirus the range of types detected in respiratory surveillance; for most a low and higher Ct value per type. For the other pathogens those detected with a Ct value below 30 during increased ILI incidence at the start of the 2018/2019 respiratory season that could not be fully explained by routinely detected viruses.

Additional were tested, CoV-HKU1 RNA from highly positive specimen, CoV-NL63 virus isolate, CoV-229E virus isolate, *Chlamydia psittaci* culture, *Chlamydia pneumoniae* culture, *Mycoplasma*

*pneumoniae* culture and the 2019 MERS-CoV QCMD EQA panel containing different dilutions of MERS-CoV in 5 specimens, CoV-OC43 (n=1) and CoV-NL63 (n=1) virus and 1 negative specimen. Except for CoV-HKU1 RNA all specimens were freshly extracted.

**Table 1.** Tests of known respiratory viruses and bacteria in clinical specimens

| Clinical specimens with known viruses/bacteria       | Number of samples tested in both assays | Ct range        |
|------------------------------------------------------|-----------------------------------------|-----------------|
| CoV-HKU1                                             | 1                                       | NA              |
| CoV-OC43                                             | 5                                       | 22.7 – 29.9     |
| CoV-229E                                             | 3                                       | 25.1 – 28.5     |
| Influenza virus A(H1N1)pdm09                         | 10                                      | 18.9 – 28.1     |
| Influenza virus A(H3N2)                              | 10                                      | 17.4 – 27.7     |
| Influenza virus B/Victoria                           | 7                                       | 19.8 – 29.9     |
| Influenza virus B/Yamagata                           | 10                                      | 23.5 – 29.1     |
| Rhinovirus                                           | 13                                      | 13.4 – 28.0     |
| Enterovirus (for types see supplementary excel file) | 18                                      | 19.6 – 28.9     |
| Respiratory syncytial virus A                        | 10                                      | 18.0 – 29.8     |
| Respiratory syncytial virus B                        | 11                                      | 21.0 – 29.4     |
| Parainfluenza virus 3                                | 1                                       | 26              |
| Human metapneumovirus                                | 3                                       | 26.7 – 27.6     |
| Adenovirus                                           | 2                                       | 21.0 – 27.3     |
| Human Bocavirus                                      | 1                                       | 35 (borderline) |
| Mycoplasma spp.                                      | 1                                       | 26.8            |
| <b>Total</b>                                         | <b>106</b>                              |                 |

#### qRT-PCR

Per clinical specimen 200 µl was extracted on a MagNA Pure 96 Instrument (Roche) using the MagNA Pure 96 DNA and Viral NA Small Volume Kit (Roche) and eluted in a volume of 50 µl. E-assay and RdRP-assay primers and probes as described in the provided protocol by Corman et al. were kindly provided by Corman and ordered from Eurogentec (Belgium). As control material SARS-CoV-1 full virus RNA (Strain Frankfurt 1, concentration 1000 copies RNA/µl) (Ct value 29-30 at 5 µl/rx) was also kindly provided by Corman.

Reaction conditions were as shown in Tables 2-4.

**Table 2.** For primers and probes provided by Corman.

|                         |       |                         |       |
|-------------------------|-------|-------------------------|-------|
| E-gene qRT PCR          | μl    | RdRP-gene qRT-PCR       | μl    |
| 4x Taqman Fast Virus MM | 5     | 4x Taqman Fast Virus MM | 5     |
| SJVE-F (10 μM)          | 1     | SJVR-F (10 μM)          | 1.5   |
| SJVE-R (10 μM)          | 1     | SJVR-R (10 μM)          | 2.0   |
| SJVE-P (10 μM)          | 0.5   | SJVR-P1/SJVR-P2 (10 μM) | 0.5   |
| PCR grade water         | 7.5   | PCR grade water         | 6     |
| Specimen nucleic acid   | 5     | Specimen nucleic acid   | 5     |
| Total volume            | 20 μl | Total volume            | 20 μl |

**Table 3.** For primers and probes obtained from Eurogentec primers and probes were premixed at a final concentration of 10 μM each except for the SJVR probes for 5 μM each.

|                         |       |                         |       |
|-------------------------|-------|-------------------------|-------|
| E-gene qRT PCR          | μl    | RdRP-gene qRT-PCR       | μl    |
| 4x Taqman Fast Virus MM | 5     | 4x Taqman Fast Virus MM | 5     |
| SJVE Mix (10 μM)        | 3     | SJVR Mix (10 μM)        | 3     |
| PCR grade water         | 7     | PCR grade water         | 7     |
| Specimen nucleic acid   | 5     | Specimen nucleic acid   | 5     |
| Total volume            | 20 μl | Total volume            | 20 μl |

**Table 4.** Amplification temperature protocol.

| PCR Program               | Segment number | Temp Target (°C) | Hold Time (sec.) | Slope (°C/sec.) | Acquisition mode | LC 480 |
|---------------------------|----------------|------------------|------------------|-----------------|------------------|--------|
| Reverse Transcription     | 1              | 50               | 900              | EXTERNAL        |                  |        |
| Denaturation/Inactivation | 1              | 95               | 120              | EXTERNAL        |                  |        |
| Denaturation              | 1              | 95               | 60               | 4.4             | None             |        |
| Amplification             | 1              | 95               | 10               | 4.4             | None             |        |
| (cycles:50)               | 2              | 58               | 30               | 2.2             | Single           |        |
| Cooling                   | 1              | 40               | 30               | 4.4             | None             |        |

Standard EAV external control qRT-PCR was performed to control for inhibitors. (IDS/VIR/F351)

**Analytic parameters:**

The detection qPCR is not used as a quantitative method. Results are scored positive/negative to confirm if a patient is infected with the novel coronavirus 2019-nCoV. Therefore all analytic parameters, used for validation of quantitative methods, are not tested.

Measurement uncertainty (Meetonzekerheid):

– not tested

Measurement trueness (Juistheid):

– tested: MERS-CoV panel and Frankfurt-1 SARS-CoV-1 control RNA (page 13)

Measurement accuracy (Accuraatheid):

– not tested, no golden standard available

Measurement precision including measurement repeatability (Precisie inclusief herhaalbaarheid):

– not tested

Analytical sensitivity (Analytische sensitiviteit):

– tested: LOD95 determination (page 12)

Analytical specificity, including interfering substances (Analytische Specificiteit met inbegrip van interfererende substanties):

– tested: clinical validation including negative controls (page 13)

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***Clinical parameters:***

Diagnostic sensitivity (Diagnostische Sensitiviteit):

– not tested, no clinical specimens of SARS-CoV-2 positive patients available

Diagnostic specificity (Diagnostische Specificiteit):

– not tested, no clinical specimens of SARS-CoV-2 positive patients available

**Results**

**LOD95 determination**

10-fold dilution series of SARS-CoV-1 control material were amplified using the reference primers and probes and the Eurogentec ordered primers and probes. Based on this dilution series a starting point for 2-fold dilution series in 4-fold was chosen, prepared and subjected to amplification to determine the LOD95.

**Table 5.** RT-qPCR results of 10-fold dilution series of SARS-CoV-1 RNA with both primer and probe sets

| Name                        | E-gene                |                      | RdRP-gene             |                      |
|-----------------------------|-----------------------|----------------------|-----------------------|----------------------|
|                             | Eurogentec (EuroE) Ct | Duitsland (DuitE) Ct | Eurogentec (EuroR) Ct | Duitsland (DuitR) Ct |
| SARS-CoV-1 undiluted        | 26.09                 | 25.15                | 27.08                 | 28.77                |
| SARS-CoV-1 10 <sup>-1</sup> | 29.18                 | 28.32                | 30.49                 | 31.48                |
| SARS-CoV-1 10 <sup>-2</sup> | 31.84                 | 30.85                | 33.23                 | 32.98                |
| SARS-CoV-1 10 <sup>-3</sup> | 33.99                 | 34.01                | 35.68                 | 34.5                 |
| SARS-CoV-1 10 <sup>-4</sup> | neg                   | 35.92                | 35.75                 | 34.97                |
| SARS-CoV-1 10 <sup>-5</sup> | neg                   | 36.12                | neg                   | neg                  |
| SARS-CoV-1 10 <sup>-6</sup> | neg                   | 35.18                | neg                   | neg                  |
| SARS-CoV-1 10 <sup>-7</sup> | neg                   | 35.26                | neg                   | neg                  |
| Neg extraction control 1    | neg                   | 34.44                | neg                   | neg                  |
| Neg extraction control 2    | neg                   | ND                   | neg                   | neg                  |

nd = not done

Eurogentec E-gene primers and probe give a cleaner picture (Figure 1). German E-gene primers and probe generate 'ghost' curves in high dilutions and negative control (Figure 1). For the RdRP-gene qRT-PCR primers and probes of both sources the qRT-PCRs generated similar results. For both qRT-PCR assays the Eurogentec primers and probes generated higher fluorescence values than the German primers and probes.


**Figure 1.** Amplification curve E-gene qRT-PCR. Blue Eurogentec primers and probe; red German primers and probe.




**Figure 2.** Amplification curve RdRP-gene qRT-PCR. Pink Eurogentec primers and probe; green German primers and probe.

Based on the 10-fold dilution series the linearity was determined and the efficiency calculated (Figure 3). Result of the LOD95% determination is shown in Figures 4 and 5.

**A**

**B**

|                | Slope  | $E = 10^{-1/\text{slope}}$ | $E (\%) = 100 \times [10^{-1/\text{slope}} - 1]$ | $R^2$  |
|----------------|--------|----------------------------|--------------------------------------------------|--------|
| Theoretical    | -3.322 | 2.00                       | 100                                              | 0.9500 |
| nCoV E-gene    | -2.382 | 2.63                       | 162.95                                           | 0.9943 |
| nCoV RdRp-gene | -2.865 | 2.23                       | 123.38                                           | 0.9979 |

**Figure 3.** Linearity (panel A) and amplification efficiency (panel B) with Eurogentec primers and probes

Linearity and efficiency were excellent.

#### LOD95% with upper and lower bound


**Figure 5:** The LOD95% copies/μL RNA of the current E-gene and RdRp-gene RT-qPCR.


**Figure 4.** 95% limit of detection curves with calculated LOD95; Eurogentec primers and probes

### Clinical specimens

Except for two of the clinical specimens and other cultured viruses, the cultured bacteria and the 2019 MERS-CoV QCMD EQA panel specimens were negative in the qRT-PCR with Eurogentec primers and probes.

One (with Ct 28.1) of 10 influenza A(H1N1)pdm09 positive specimens and one (with Ct 29.3) of 11 RSV-B positive specimens were positive in the E-gene qRT-PCR only with high Ct values of 37.7 and 39.1 respectively (Figure 5). Both specimens have been repeated in the E-gene RT-PCR with annealing temperatures of 58°C and 60°C and turned out negative. None of the other A(H1N1)pdm09 and RSV-B positive specimens with higher viral load generated a positive signal in the E-gene qRT-PCR.


Figure 6. Amplification curves of clinical specimens with E-gene qRT-PCR

### Discussion

Except for two of the clinical specimens and other cultured viruses, the cultured bacteria and the 2019 MERS-CoV QCMD EQA panel specimens were negative in the qRT-PCR with Eurogentec primers and probes. Both specimens have been repeated in the E-gene RT- and turned out negative. The measurement trueness is 100%.

Based on the 10-fold dilution series in both E-gene and RdRP-Gene the efficiency excellent.

Based on the 2-fold dilution series in both E-gene and RdRP-Gene the analytical sensitivity is below 1 copy per  $\mu$ l input RNA and 3.2 and 3.7 copies RNA/reaction, respectively.

### Conclusion

Based on the results of the measurement trueness, efficiency and analytical specificity and analytical sensitivity the test can be implemented to detect the new coronavirus 2019-nCoV.

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