

1 Supplementary tables

2

Table 1. Sample panel used to validate the performance of the ELISA/Liaison and the antibody RDT assays for SARS-CoV-2				
Sensitivity				
Country	Sample source	Infection	No. samples Liaison/ELISA (RDT*)	Post symptom onset range
Netherlands	RT-PCR confirmed SARS-CoV-2	Mild	15/15 (15)	1-24
		Severe	78/78 (78)	1-43
Specificity				
Netherlands	Healthy blood donors	NA	20/20 (11)	NA
Netherlands	Non-CoV respiratory infections*	Adeno virus	6/6 (1)	2-4weeks
		HMPV	9/9 (3)	2-4weeks
		Flu A	13/13 (4)	2-4weeks
		Flu B	6/6 (4)	2-4weeks
		RSV A	5/5 (4)	2-4weeks
		RSV B	4/4 (4)	2-4weeks
		CMV	5/5 (2)	2-4weeks
		EBV	7/7 (3)	2-4weeks
		<i>M.pneumoniae</i>	1/1 (1)	2-4weeks
		Rhinovirus	9/9 (2)	2-4weeks
		Bocavirus	2/2 (0)	2-4weeks
		Parainfluenza 1/3	8/8 (0)	2-4weeks
		Enterovirus	2/2 (0)	2-4weeks
Netherlands	hCoV infections	HCoV 229E	19/19 (6)	2-4weeks
		HCoV-NL63	17/17 (7)	2-4weeks
		HCoV-OC43	39/39 (9)	2-4weeks
Soudan/Qatar		MERS	8/8 (3)	2-4weeks

3

4 * numbers were limited due to RDT kit availability

5

6

7

8

9

10

11

12

13

14

15 **Table 2. A summary of the performance characteristics of the Wantai Ig and IgM ELISA, Euroimmun**
 16 **IgG and IgA ELISA and DiaSorin Liason platform.**

17

		Wantai Ig	Wantai IgM	Euroimmun IgG*	Euroimmun IgA*	Liaison	
Specificity	n/N	149/150	148/150	160/161	151/161	68/69	
	Value	0.9933	0.9867	0.9938	0.9379	0.9855	
	(95% CI)	{0.9632 - 0.9997}	{0.9527 - 0.9976}	{0.9657 - 0.9997}	{0.8895 - 0.9659}	{0.9224 - 0.9993}	
Overall	Sensitivity	n/N	75/76	68/76	62/76	74/76	39/53
		Value	0.9868	0.8947	0.8158	0.9737	0.7358
		(95% CI)	{0.9292 - 0.9993}	{0.8058 - 0.9457}	{0.7142 - 0.8870}	{0.9090 - 0.9953}	{0.6042 - 0.8356}
Positive Predictive Value	Value	0.9868	0.9714	0.9841	0.881	0.975	
	(95% CI)	{0.9292 - 0.9993}	{0.9017 - 0.9949}	{0.9154 - 0.9992}	{0.7946 - 0.9340}	{0.8712 - 0.9987}	
Negative Predictive Value	Value	0.9933	0.9487	0.9195	0.9869	0.8293	
	(95% CI)	{0.9632 - 0.9997}	{0.9021 - 0.9738}	{0.8695 - 0.9515}	{0.9536 - 0.9977}	{0.7336 - 0.8955}	
Likelihood Ratio		148	67.11	131.3	15.68	50.77	
>14 dps	Sensitivity	n/N	26/26	19/26	25/26	26/26	17/18
		Value	1	0.7308	0.9615	1	0.9444
		(95% CI)	{0.8713 - 1.000}	{0.5392 - 0.8630}	{0.8111 - 0.9980}	{0.8713 - 1.000}	{0.7424 - 0.9972}
Positive Predictive Value	Value	0.963	0.9048	0.9615	0.7222	0.9444	
	(95% CI)	{0.8172 - 0.9981}	{0.7109 - 0.9831}	{0.8111 - 0.9980}	{0.5601 - 0.8415}	{0.7424 - 0.9972}	
Negative Predictive Value	Value	1	0.9548	0.9938	1	0.9855	
	(95% CI)	{0.9749 - 1.000}	{0.9097 - 0.9780}	{0.9657 - 0.9997}	{0.9752 - 1.000}	{0.9224 - 0.9993}	
Likelihood Ratio		150	54.81	154.8	16.1	65.17	

18

19

20

21

22

23

24

25

26

27

28

29

30

Table 3. Sensitivity and specificity of the three tested RDTs (CI 95%) compared to PRNT50 positivity: whole sample set and samples >14 days post symptom onset (borderline PRNT50 values counted as positive).						
	Cellex (n=93)		InTec (n=93)		Orient gene (n=90)	
	IgM	IgG	IgM	IgG	IgM	IgG
Compared to PRNT50 positivity						
Sensitivity-overall	87.36% (78.50% to 93.52%)	84.44% (75.28% to 91.23%)	88.37% (79.65% to 94.28%)	95.00% (87.69% to 98.62%)	89.41% (80.85% to 95.04%)	91.57% (83.39% to 96.54%)
Sensitivity >14 DPO (n=25)	96.15% (80.36% to 99.90%)	96.15% (80.36% to 99.90%)	88.37% (79.65% to 94.28%)	96.15% (80.36% to 99.90%)	80.65% (62.53% to 92.55%)	100.00% (86.28% to 100.00%)
	IgM and IgG		IgM and IgG		IgM and IgG	
	96.15% (80.36% to 99.90%)		96.15% (80.36% to 99.90%)		100.00% (86.28% to 100.00%)	
Specificity - overall*	80.95% (58.09% to 94.55%)	85.00% (62.11% to 96.79%)	73.91% (51.59% to 89.77%)	77.27% (54.63% to 92.18%)	100.00% (80.49% to 100.00%)	100.00% (80.49% to 100.00%)

31 * number of sera tested for specificity based on availability of test kit; n=44, 64 and 9 respective of

32 Cellex, inTec and Orient Gene

33

34

35

36

37

38

39

40
41