

Corona App from Münster

Speed-up testing &
Improve patient follow-up

5.1.2e

5.1.2e

@hippokrates-it.de



Workflow of Corona App: Speed-up of testing

- Person is eligible for COVID-19 test
- Download of app to smartphone
- Standardized questionnaire is filled in before testing
- 2D-Barcode is generated on Smartphone
- Barcode is scanned into IT system of test center

8:22

My Patient Data

APR 10, 2020 8:22 AM

CORONA MED. HISTORY

This app is used to prepare Your medical encounter regarding a possible corona infection. Please answer the questions carefully in advance and then speak to your doctor. The app generates a barcode from your data, which is read into your doctor's computer. NO data is transferred to the Internet. This helps to speed-up the diagnosis process.

Symptoms

Do you have a cough? (*)
 yes no

8:22

Symptoms

Do you have a cough? (*)
 yes no

Have you got the sniffles? (*)
 yes no

Do you have hoarseness? (*)
 yes no

Do you have a fever or increased body temperature? (*)
 yes no

What body temperature have you currently measured?
 39 Degree Celsius

Do You have a disturbed sense of smell? (*)
 yes no

Do you have a disturbed sense of taste? (*)
 yes no

Do You have diarrhea? (*)

8:23

Do you have a disturbed sense of taste? (*)
 yes no

Do You have diarrhea? (*)
 yes no

Contact

Did you have contact to people with corona infection in the past 14 days? (*)
 yes no

Contact person
 Mr ABC

Do you live/reside/work in a health care facility (nursing home/hospital) where respiratory diseases have frequently occurred?
 yes no

Risk area

Have you been in a risk area for corona infection during the past 14 days or have you had contact with people who have been

8:24

Risk area

Have you been in a risk area for corona infection during the past 14 days or have you had contact with people who have been there?

Risk areas are: Egypt, France, Iran, Italy, Netherlands, Austria, Switzerland, Spain, South Korea: Daegue and Gyeongsangbuk-do Province, USA, Great Britain and Northern Ireland (see ecdc.europa.eu for the current list of risk areas) (*)

yes no

Where were you or your contacts?

Netherlands

General information

Do you have other current diseases? (*)

yes no

Which current diseases?

MDS

8:25

Birth date (*)

1.1.1950

Phone number (*)

1234567

Do you agree to be contacted by the doctor/clinic at a later date?

yes no

Additional information

Study Patient

BACK **STORE**

8:25

+ **BARCODE**

10-Apr-20 08:22 Corona Med. History :


First name:Werner Last name:Test Birth date:1.1.1950 Phone:1234567 Symptoms: cough, sniff, hoarseness, fever, temperature: 39 degree Celsius, dysosmia, dysgeusia, diarrhoa. Contact to infected people: YES. Contact person: Mr ABC. Contact to respiratory diseases in healthcare facilities: YES. Stay in risk area: YES. location: Netherlands. General information: Current diseases: MDS. Chronic diseases: Diabetes. Medication: Insulin. Study Patient FU yes



My Patient Data



Entry 1-1



GO BACK

SCANNER



Workflow of Corona App: Improve patient follow-up

- Person is sent into quarantine
- App sends daily reminders to document vital signs & symptoms
- 2D-barcode is generated on smartphone
- Barcode is scanned into IT system of GP or hospital
- Physician can review disease course

8:26

FOLLOW-UP

Please document the course of your illness regularly. If you feel acutely worse, please contact your doctor and report those data.

Vital signs

What's your temperature?
 39 , 0 Degree Celsius (*)

Pulse /min (*)

Please enter your systolic blood pressure (upper value of blood pressure measurement)
 mmHg (*)

Symptoms

On a scale from 0 to 10, how severe is your cough? (*)

0 No cough

8:26

Symptoms

On a scale from 0 to 10, how severe is your cough? (*)

0 No cough

1

2

3

4

5

6

7

8

9

10 Extreme cough

Are you experiencing shortness of breath? (*)

Not troubled by breathless except on strenuous exercise

Short of breath when hurrying on a level or when walking up a slight hill

Walks slower than most people on the level, stops after a mile or so, or stops after 15 minutes walking at own pace

Stops for breath after walking 100 yards, or after a few minutes on level ground

Too breathless to leave the house, or

8:27

If you have a pulse oximeter: What's the oxygen saturation (SpO2)?

Have you got the sniffles?
 yes no

Do you have diarrhoea?
 yes no

Do you have a disturbed sense of smell?
 yes no

Do you have a disturbed sense of taste?
 yes no

If you have other symptoms, what are they?
Headache

Medication

If you take medication regularly: Have you taken your medication today?
 yes no

If you took any extra medicine today: What

Corona App: Key characteristics

- Scalability: No additional software installation & compatible with ANY clinical software
- GDPR-compliant: no data transfer via Internet
- No Geo-Tracking
- Minimal setup cost: only 1 standard Barcode scanner per test center / GP / hospital
- App is free of charge, no advertisement

Corona App: Value proposition

- Speed-up of testing up to 20%: less time for medical history taking
- Contact time is reduced
- Patient triage is supported by regular documentation
- Tested in routine operation for 3 weeks
- Available in 8 languages: German, English, Spanish, Portuguese, Italian, Dutch, French, Turkish
- Ready for use now

Corona App from Münster

Speed-up testing &
Improve patient follow-up

More infos:

<https://hippokrates-it.de/corona/>

