



85.0
175.0



**恒安集团**  
HENGAN



Disposable Surgical Mask

**How to use**

1. Clean your hands
2. Check the medical mask for defects
3. Check production and expiration date
4. Orient the top of the mask properly
5. Ensure the proper side of the mask faces outwards
6. Place the mask on your face
7. Adjust the nose piece
8. Fit the mask to your face and under your chin

For clinical staff to wear during invasive operations.  
Wear to provide protection for patients treated by medical staff.  
Stops the spread of blood and body fluids.

**Specifications:**

- 17,5 cm x 9,5 cm (+/- 0,5 cm)
- Earloops
- Adjustable nose piece
- Disposable / single use only

10PCS

**Shelf life:**  
18 months from date of manufacture

**Storage conditions:**  
Keep in a dry environment and away from excessive heat

**Manufactured by:**  
Jinjiang Hengan Household Tissue Products Co., Ltd.  
Wuli Industrial Zone, Anhai Town, Jinjiang City, Fujian, China  
www.hengan.com

**Imported by:**  
Mediq Nederland B.V.  
On behalf of Landelijk Consortium Hulpmiddelen  
Rijnzathe 19, 3454 PV Utrecht, The Netherlands


BFE ≥ 95 % / PFE ≥ 30 %

【医疗器械注册证编号】：  
闽械注准20202140015 (临时)

【生产许可证编号】：  
闽药监械生产许20200437号 (临时)

**Standard:** YY0469-2011

**Model Number:** XCJ010



6 903244 681780